

CD3EAP Rabbit pAb

Catalog No.: A16099 **2 Publications**

Basic Information

Observed MW

70kDa

Calculated MW

55kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse

Background

Enables RNA binding activity. Predicted to be involved in transmembrane receptor protein tyrosine kinase signaling pathway. Predicted to act upstream of or within rRNA transcription. Located in cytosol; mitochondrion; and nuclear lumen.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

10849

Swiss Prot

O15446

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human CD3EAP (NP_036231.1).

Synonyms

ASE1; CAST; ASE-1; PAF49; RPA34; CD3EAP

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using CD3EAP Rabbit pAb (A16099) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 1s.