

TMC5 Rabbit pAb

Catalog No.: A16145

Basic Information

Observed MW

114kDa

Calculated MW

115kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Mouse, Rat

Background

Predicted to enable mechanosensitive ion channel activity. Predicted to be involved in ion transmembrane transport. Located in extracellular exosome.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

79838

Swiss Prot

Q6UXY8

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 676-760 of human TMC5 (NP_079056.2).

Synonyms

TMC5

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using TMC5 Rabbit pAb (A16145) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.