

BAZ1A Rabbit pAb

Catalog No.: A16495

Basic Information

Observed MW

Refer to figures

Calculated MW

179kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

The BAZ1A gene encodes the accessory subunit of the ATP-dependent chromatin assembly factor (ACF), a member of the ISWI ('imitation switch') family of chromatin remodeling complexes (summarized by Racki et al., 2009 [PubMed 20033039]).

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

11177

Swiss Prot

Q9NRL2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1347-1556 of human BAZ1A (NP_038476.2).

Synonyms

ACF1; WALp1; hACF1; WCRF180; BAZ1A

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.