

# FASTKD2 Rabbit pAb

Catalog No.: A16499

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

81kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human

## Background

This gene encodes a protein that is localized in the mitochondrial inner compartment and that may play a role in mitochondrial apoptosis. Nonsense mutations have been reported to result in cytochrome c oxidase deficiency.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

22868

### Swiss Prot

Q9NYY8

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 439-710 of human FASTKD2 (NP\_055744.2).

### Synonyms

COXPD44; KIAA0971; FASTKD2

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.