

MON2 Rabbit pAb

Catalog No.: A16501

Basic Information

Observed MW

Refer to figures

Calculated MW

190kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Predicted to be involved in Golgi to endosome transport. Located in early endosome membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

23041

Swiss Prot

Q7Z3U7

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human MON2 (NP_001265399.1).

Synonyms

MON2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.