

# ENDOD1 Rabbit pAb

Catalog No.: A16502 **1 Publications**

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

55kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human, Mouse, Rat

## Background

Predicted to enable endonuclease activity; metal ion binding activity; and nucleic acid binding activity. Predicted to be involved in nucleic acid phosphodiester bond hydrolysis. Located in extracellular exosome and membrane.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

23052

### Swiss Prot

O94919

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 22-330 of human ENDOD1 (NP\_055851.1).

### Synonyms

ENDOD1

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.