

# PTPLAD1 Rabbit pAb

Catalog No.: A16531

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

43kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse

## Background

Enables 3-hydroxyacyl-CoA dehydratase activity and enzyme binding activity. Involved in fatty acid biosynthetic process; positive regulation by virus of viral protein levels in host cell; and positive regulation of viral genome replication. Located in endoplasmic reticulum and nuclear membrane.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

51495

### Swiss Prot

Q9P035

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human PTPLAD1 (NP\_057479.2).

### Synonyms

BIND1; B-IND1; HSPC121; PTPLAD1

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.