

# TDRD9 Rabbit pAb

Catalog No.: A16592

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

156kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse, Rat

## Background

Predicted to enable RNA binding activity. Involved in spermatogenesis. Located in cytoplasm and nucleus. Implicated in spermatogenic failure 30.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

122402

### Swiss Prot

Q8NDG6

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 300-600 of human TDRD9 (NP\_694591.2).

### Synonyms

HLS; SPNE; HIG-1; NET54; SPGF30; C14orf75; TDRD9

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.