

SYNPR Rabbit pAb

Catalog No.: A16595

Basic Information

Observed MW

Refer to figures

Calculated MW

29kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

Background

Predicted to be located in neuron projection and synaptic vesicle. Predicted to be integral component of membrane. Predicted to be active in synaptic vesicle membrane. Predicted to be integral component of synaptic vesicle membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

132204

Swiss Prot

Q8TBG9

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

SPO; SYNPR

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.