

RNF152 Rabbit pAb

Catalog No.: A16609

Basic Information

Observed MW

Refer to figures

Calculated MW

22kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Rat

Background

Enables small GTPase binding activity and ubiquitin protein ligase activity. Involved in several processes, including cellular response to amino acid starvation; negative regulation of TORC1 signaling; and protein polyubiquitination. Located in lysosomal membrane. Is integral component of organelle membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

220441

Swiss Prot

Q8N8N0

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RNF152 (NP_775828.1).

Synonyms

RNF152

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.