

NFIA Rabbit pAb

Catalog No.: A16719

Basic Information

Observed MW

Refer to figures

Calculated MW

56kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Rat

Background

This gene encodes a member of the NF1 (nuclear factor 1) family of transcription factors. Multiple transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:200 - 1:2000

Immunogen Information

Gene ID

4774

Swiss Prot

Q12857

Immunogen

A synthetic Peptide of human NFIA .

Synonyms

CTF; NF1-A; NFI-A; NFI-L; BRMUTD; NF-I/A; DEL1P32P31; C1DELp32p31; NFIA

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.