

APOD Rabbit pAb

Catalog No.: A16758

Basic Information

Observed MW

Refer to figures

Calculated MW

21kDa

Category

Primary antibody

Applications

ELISA, IF/ICC

Cross-Reactivity

Human

Background

This gene encodes a component of high density lipoprotein that has no marked similarity to other apolipoprotein sequences. It has a high degree of homology to plasma retinol-binding protein and other members of the alpha 2 microglobulin protein superfamily of carrier proteins, also known as lipocalins. This glycoprotein is closely associated with the enzyme lecithin:cholesterol acyltransferase - an enzyme involved in lipoprotein metabolism.

Recommended Dilutions

IF/ICC 1:50 - 1:100

Immunogen Information

Gene ID

347

Swiss Prot

P05090

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 21-189 of human APOD (NP_001638.1).

Synonyms

APOD

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.