

TRAPPC8 Rabbit pAb

Catalog No.: A17099

Basic Information

Observed MW

Refer to figures

Calculated MW

161kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Rat

Background

Involved in Golgi organization. Part of TRAPP complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

22878

Swiss Prot

Q9Y2L5

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

GSG1; TRS85; HsT2706; KIAA1012; TRAPPC8

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.