

RAPGEF6 Rabbit pAb

Catalog No.: A17162

Basic Information

Observed MW

Refer to figures

Calculated MW

179kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables several functions, including GTP-dependent protein binding activity; guanyl-nucleotide exchange factor activity; and phosphatidic acid binding activity. Involved in microvillus assembly; positive regulation of GTPase activity; and protein localization to plasma membrane. Located in several cellular components, including apical plasma membrane; centrosome; and endocytic vesicle.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

51735

Swiss Prot

Q8TEU7

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

RAGEF2; PDZGEF2; KIA001LB; PDZ-GEF2; RA-GEF-2; RAPGEF6

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.