

TCEAL6 Rabbit pAb

Catalog No.: A17266

Basic Information

Observed MW

Refer to figures

Calculated MW

22kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Predicted to enable WW domain binding activity. Predicted to be active in nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

158931

Swiss Prot

Q6IPX3

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

WEX2; Tceal3; TCEAL6

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.