

Kctd19 Rabbit pAb

Catalog No.: A17274 **2 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

105kDa

Category

Primary antibody

Applications

ELISA,IF/ICC

Cross-Reactivity

Mouse, Rat

Background

Predicted to be involved in protein homooligomerization. Orthologous to human KCTD19 (potassium channel tetramerization domain containing 19).

Recommended Dilutions

IF/ICC 1:50 - 1:200

Immunogen Information

Gene ID

279499

Swiss Prot

Q562E2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of mouse Kctd19 (NP_001288102.1).

Synonyms

4922504H04Rik; Kctd19

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.