

C10orf62 Rabbit pAb

Catalog No.: A17284

Basic Information

Observed MW

Refer to figures

Calculated MW

25kDa

Category

Primary antibody

Applications

ELISA, IHC-P

Cross-Reactivity

Mouse, Rat

Background

Recommended Dilutions

IHC-P 1:50 - 1:200

Immunogen Information

Gene ID

414157

Swiss Prot

Q5T681

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human C10orf62 (NP_001009997.2).

Synonyms

bA548K23.1; C10orf62

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Immunohistochemistry analysis of C10orf62 in paraffin-embedded rat testis using C10orf62 Rabbit pAb (A17284) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry analysis of C10orf62 in paraffin-embedded mouse testis using C10orf62 Rabbit pAb (A17284) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.