

MIOX Rabbit pAb

Catalog No.: A17287

Basic Information

Observed MW

Refer to figures

Calculated MW

46kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

Background

Predicted to enable aspartic-type peptidase activity; glutamic-type peptidase activity; and metalloendopeptidase activity. Predicted to be involved in several processes, including negative regulation of binding activity of sperm to zona pellucida; positive regulation of protein processing; and prevention of polyspermy. Predicted to be located in cortical granule and plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

431705

Swiss Prot

Q6HA08

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

SAS1B; OOMD11; OZEMA11; MIOX

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.