

SIDT2 Rabbit pAb

Catalog No.: A17317

Basic Information

Observed MW

Refer to figures

Calculated MW

94kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Rat

Background

Predicted to enable several functions, including AP-1 adaptor complex binding activity; AP-2 adaptor complex binding activity; and RNA transmembrane transporter activity. Involved in RNA transport. Located in lysosomal membrane and plasma membrane.

Recommended Dilutions

Immunogen Information

Gene ID

51092

Swiss Prot

Q8NBJ9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 19-290 of human SIDT2 (NP_001035545.1).

Synonyms

CGI-40; SIDT2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.