

UBAP2L Rabbit pAb

Catalog No.: A17371

Basic Information

Observed MW

115kDa

Calculated MW

115kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Enables RNA binding activity. Involved in binding activity of sperm to zona pellucida and stress granule assembly. Acts upstream of or within hematopoietic stem cell homeostasis. Part of PcG protein complex.

Recommended Dilutions

WB 1:1000 - 1:5000

Immunogen Information

Gene ID

9898

Swiss Prot

Q14157

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 868-1087 of human UBAP2L (NP_055662.3).

Synonyms

NICE4; NICE-4; UBAP2L

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from A-549 cells, using UBAP2L Rabbit pAb (A17371) at 1:2000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 180s.