

ZNHIT3 Rabbit pAb

Catalog No.: A17375

Basic Information

Observed MW

Refer to figures

Calculated MW

18kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Predicted to enable thyroid hormone receptor binding activity. Predicted to be involved in box C/D snoRNP assembly; maturation of LSU-rRNA from tricistronic rRNA transcript (SSU-rRNA, 5.8S rRNA, LSU-rRNA); and snoRNA localization. Located in cytoplasm and nucleus. Implicated in PEHO syndrome.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

9326

Swiss Prot

Q15649

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 35-100 of human ZNHIT3 (NP_004764.1).

Synonyms

Hit1; PEHO; TRIP3; ZNHIT3

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.