

TAOK1 Rabbit pAb

Catalog No.: A17405

Basic Information

Observed MW

Refer to figures

Calculated MW

116kDa

Category

Primary antibody

Applications

WB,IF/ICC,IP

Cross-Reactivity

Human, Mouse

Background

Enables alpha-tubulin binding activity; beta-tubulin binding activity; and kinase activity. Involved in several processes, including mitotic G2 DNA damage checkpoint signaling; negative regulation of microtubule depolymerization; and positive regulation of JNK cascade. Located in microtubule cytoskeleton and perinuclear region of cytoplasm.

Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:200
IP	1:50 - 1:200

Immunogen Information

Gene ID

57551

Swiss Prot

Q7L7X3

Immunogen

Recombinant Protein of human TAOK1.

Synonyms

DDIB; PSK2; TAO1; KFC-B; MARKK; PSK-2; hKFC-B; hTAOK1; MAP3K16; TAOK1

Contact

☎	400-999-6126
✉	cn.market@abclonal.com.cn
🌐	www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.