

WWC2 Rabbit pAb

Catalog No.: A17410

Basic Information

Observed MW

Refer to figures

Calculated MW

134kDa

Category

Primary antibody

Applications

WB, IF/ICC, IP

Cross-Reactivity

Human, Mouse

Background

This gene encodes a member of the WW-and-C2-domain-containing family of proteins. Members of this family have two N-terminal WW domains that mediate binding to target proteins harboring L/PPxY motifs, an internal C2 domain for membrane association, and C-terminal alpha protein kinase C binding sites and class III PDZ domain-interaction motifs. Proteins of this family are able to form homo- and heterodimers and to modulate hippo pathway signaling.

Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:200
IP	1:50 - 1:200

Immunogen Information

Gene ID

80014

Swiss Prot

Q6AWC2

Immunogen

Recombinant Protein of human WWC2.

Synonyms

BOMB; WWC2

Contact

☎	400-999-6126
✉	cn.market@abclonal.com.cn
🌐	www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.