

# WWC3 Rabbit pAb

Catalog No.: A17411

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

123kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human

## Background

This gene encodes a member of the WWC family of proteins, which also includes the WWC1 (KIBRA) gene product and the WWC2 gene product. The protein encoded by this gene includes a C2 domain, which is known to mediate homodimerization in the related WWC1 gene product.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

55841

### Swiss Prot

Q9ULE0

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 580-760 of human WWC3 (NP\_056506.2).

### Synonyms

BM042; WWC3

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.