

EPG5 Rabbit pAb

Catalog No.: A17456

Basic Information

Observed MW

Refer to figures

Calculated MW

292kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

This gene encodes a large coiled coil domain-containing protein that functions in autophagy during starvation conditions. Mutations in this gene cause Vici syndrome.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

57724

Swiss Prot

Q9HCE0

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 819-947 of human EPG5 (NP_066015.2).

Synonyms

HEEW1; VICIS; KIAA1632; EPG5

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.