

GK2 Rabbit pAb

Catalog No.: A17487

Basic Information

Observed MW

Refer to figures

Calculated MW

61kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Rat

Background

Predicted to enable glycerol kinase activity. Predicted to be involved in several processes, including glycerol metabolic process; glycerol-3-phosphate biosynthetic process; and triglyceride metabolic process. Located in extracellular exosome.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

2712

Swiss Prot

Q14410

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

GKP2; GKTA; GK2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.