

SH3BGRL Rabbit pAb

Catalog No.: A17530

Basic Information

Observed MW

16kDa

Calculated MW

13kDa

Category

Primary antibody

Applications

WB,IF/ICC,ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to enable SH3 domain binding activity. Located in extracellular exosome.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

6451

Swiss Prot

O75368

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

SH3BGR; HEL-S-115; SH3BGRL

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG



Purification

Affinity purification

Storage

Store at -20°C

Validation Data



Immunofluorescence analysis of U-2 OS cells using SH3BGRL Rabbit pAb (A17530) at dilution of 1:100. Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.