

CRLF1 Rabbit pAb

Catalog No.: A17589 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

46kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Mouse, Rat

Background

This gene encodes a member of the cytokine type I receptor family. The protein forms a secreted complex with cardiotrophin-like cytokine factor 1 and acts on cells expressing ciliary neurotrophic factor receptors. The complex can promote survival of neuronal cells. Mutations in this gene result in Crisponi syndrome and cold-induced sweating syndrome.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

9244

Swiss Prot

O75462

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 135-350 of human CRLF1 (NP_004741.1).

Synonyms

CLF; NR6; CISS; CISS1; CLF-1; zcytor5; CRLF1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.