

RHOBTB3 Rabbit pAb

Catalog No.: A17638 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

69kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse

Background

Enables ATP binding activity and small GTPase binding activity. Involved in retrograde transport, endosome to Golgi. Located in extracellular exosome.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

22836

Swiss Prot

O94955

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human RHOBTB3 (NP_055714.3).

Synonyms

RHOBTB3

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.