

# ANKRD6 Rabbit pAb

Catalog No.: A17642

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

80kDa

### Category

Primary antibody

### Applications

WB, ELISA

### Cross-Reactivity

Rat

## Background

Predicted to be involved in negative regulation of canonical Wnt signaling pathway and positive regulation of JNK cascade. Predicted to act upstream of or within positive regulation of Wnt signaling pathway, planar cell polarity pathway. Located in intracellular membrane-bounded organelle.

## Recommended Dilutions

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

22881

### Swiss Prot

Q9Y2G4

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

ANKRD6

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.