

WDTC1 Rabbit pAb

Catalog No.: A17650

Basic Information

Observed MW

Refer to figures

Calculated MW

76kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to enable enzyme inhibitor activity; histone binding activity; and histone deacetylase binding activity. Predicted to be involved in negative regulation of fatty acid biosynthetic process. Predicted to act upstream of or within several processes, including cellular response to insulin stimulus; glucose metabolic process; and negative regulation of transcription by RNA polymerase II. Predicted to be located in nucleus. Predicted to be part of Cul4-RING E3 ubiquitin ligase complex. Predicted to be active in cytoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

23038

Swiss Prot

Q8N5D0

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ADP; DCAF9; WDTC1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.