

# KDM4C Rabbit pAb

Catalog No.: A17652

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

120kDa

### Category

Primary antibody

### Applications

ELISA

### Cross-Reactivity

Human, Rat

## Background

This gene is a member of the Jumonji domain 2 (JMJD2) family. The encoded protein is a trimethylation-specific demethylase, and converts specific trimethylated histone residues to the dimethylated form. This enzymatic action regulates gene expression and chromosome segregation. Chromosomal aberrations and changes in expression of this gene may be found in tumor cells. Alternative splicing results in multiple transcript variants.

## Recommended Dilutions

## Immunogen Information

### Gene ID

23081

### Swiss Prot

Q9H3R0

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 350-570 of human KDM4C (NP\_055876.2).

### Synonyms

GASC1; JHDM3C; JMJD2C; TDRD14C; KDM4C

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.