

COBL Rabbit pAb

Catalog No.: A17661

Basic Information

Observed MW

Refer to figures

Calculated MW

136kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Rat

Background

This gene encodes a protein that contains WH2 domains (WASP, Wiskott-Aldrich syndrome protein, homology domain-2) that interact with actin. The encoded actin regulator protein is required for growth and assembly of brush border microvilli that play a role in maintaining intestinal homeostasis. A similar protein in mouse functions in midbrain neural tube closure. A pseudogene of this gene is located on chromosome X.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

23242

Swiss Prot

O75128

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

COBL

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.