

HACL1 Rabbit pAb

Catalog No.: A17674

Basic Information

Observed MW

Refer to figures

Calculated MW

64kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse

Background

Enables several functions, including 2-hydroxy-3-methylhexadecanoyl-CoA lyase activity; ATP binding activity; and cation binding activity. Involved in fatty acid alpha-oxidation; phytanic acid metabolic process; and protein targeting to peroxisome. Located in nucleoplasm and peroxisome.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

26061

Swiss Prot

Q9UJ83

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

HPCL; HPCL2; PHYH2; 2-HPCL; HACL1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.