

FAM127B Rabbit pAb

Catalog No.: A17676

Basic Information

Observed MW

18kDa

Calculated MW

13kDa

Category

Primary antibody

Applications

ELISA, WB, IF/ICC

Cross-Reactivity

Human, Mouse

Background

Predicted to be active in nucleolus.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:50 - 1:100

Immunogen Information

Gene ID

26071

Swiss Prot

Q9BWD3

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human FAM127B (NP_001071640.1).

Synonyms

CXX1b; MAR8A; SIRH6; FAM127B

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification



Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Immunofluorescence analysis of L929 cells using FAM127B Rabbit pAb (A17676) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.