

MBIP Rabbit pAb

Catalog No.: A17707

Basic Information

Observed MW

39kDa

Calculated MW

39kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Rat

Background

Enables identical protein binding activity and protein kinase inhibitor activity. Involved in histone H3 acetylation; positive regulation of JNK cascade; and positive regulation of gene expression. Located in cytosol; nucleolus; and nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

51562

Swiss Prot

Q9NS73

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-344 of human MBIP (NP_057670.2).

Synonyms

MBIP

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from Rat brain, using MBIP Rabbit pAb (A17707) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 15s.