

TBC1D16 Rabbit pAb

Catalog No.: A17821

Basic Information

Observed MW

Refer to figures

Calculated MW

86kDa

Category

Primary antibody

Applications

IF/ICC,ELISA

Cross-Reactivity

Mouse, Rat

Background

Enables GTPase activator activity. Involved in regulation of receptor recycling. Located in cytosol and early endosome.

Recommended Dilutions

IF/ICC 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

125058

Swiss Prot

Q8TBPO

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

TBC1D16

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.