

ADAD2 Rabbit pAb

Catalog No.: A17840

Basic Information

Observed MW

Refer to figures

Calculated MW

62kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable double-stranded RNA adenosine deaminase activity; double-stranded RNA binding activity; and tRNA-specific adenosine deaminase activity. Predicted to be involved in RNA processing and adenosine to inosine editing. Predicted to be active in cytoplasm and nucleolus.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

161931

Swiss Prot

Q8NCV1

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

TENRL; ADAD2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.