

# KSR2 Rabbit pAb

Catalog No.: A17866

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

108kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human, Rat

## Background

Predicted to enable MAP-kinase scaffold activity; mitogen-activated protein kinase kinase binding activity; and protein kinase activity. Predicted to be involved in Ras protein signal transduction; calcium-mediated signaling; and positive regulation of cold-induced thermogenesis. Predicted to act upstream of or within positive regulation of MAPK cascade. Predicted to be active in cytosol and plasma membrane.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

283455

### Swiss Prot

Q6VAB6

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 780-890 of human KSR2 (NP\_775869.3).

### Synonyms

KSR2

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.