

# PJA2 Rabbit pAb

Catalog No.: A17871

## Basic Information

### Observed MW

Refer to Figures

### Calculated MW

78kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

Enables protein kinase A catalytic subunit binding activity; protein kinase A regulatory subunit binding activity; and ubiquitin-protein transferase activity. Involved in several processes, including protein ubiquitination; regulation of macrophage activation; and regulation of signal transduction. Located in cytoplasm; intermediate filament cytoskeleton; and plasma membrane.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

9867

### Swiss Prot

O43164

### Immunogen

A synthetic peptide of human PJA2.

### Synonyms

RNF131; Neurodap1; PJA2

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.