

PSBA Rabbit pAb

Catalog No.: A17956 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Arabidopsis thaliana

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

844802

Swiss Prot

Immunogen

A synthetic peptide of Arabidopsis thaliana PSBA.

Synonyms

PSBA

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.