

FGF2 Rabbit pAb

Catalog No.: A17989 **1 Publications**

Basic Information

Observed MW

18kDa,22kDa,24kDa/24kDa

Calculated MW

31kDa

Category

Primary antibody

Applications

ELISA,WB

Cross-Reactivity

Mouse

Background

The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members bind heparin and possess broad mitogenic and angiogenic activities. This protein has been implicated in diverse biological processes, such as limb and nervous system development, wound healing, and tumor growth. The mRNA for this gene contains multiple polyadenylation sites, and is alternatively translated from non-AUG (CUG) and AUG initiation codons, resulting in five different isoforms with distinct properties. The CUG-initiated isoforms are localized in the nucleus and are responsible for the intracrine effect, whereas, the AUG-initiated form is mostly cytosolic and is responsible for the paracrine and autocrine effects of this FGF.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

2247

Swiss Prot

P09038

Immunogen

Recombinant Protein corresponding to a sequence within amino acids 143-288 of FGF2 (NP_001997.5).

Synonyms

BFGF; FGFB; FGF-2; HBGF-2; FGF2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of lysates from Mouse lung using FGF2 Rabbit pAb(A17989) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 10s.