

C10orf67 Rabbit pAb

Catalog No.: A17993

Basic Information

Observed MW

64kDa

Calculated MW

64kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to be located in mitochondrion.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

256815

Swiss Prot

Q8IYJ2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 68-243 of C10orf67 (NP_001338235.1).

Synonyms

C10orf115; LINC01552; C10orf67

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from Rat testis, using C10orf67 Rabbit pAb (A17993) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 180s.