

STOX1 Rabbit pAb

Catalog No.: A18137

Basic Information

Observed MW

Refer to figures

Calculated MW

111kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Enables RNA polymerase II transcription regulatory region sequence-specific DNA binding activity. Involved in several processes, including positive regulation of G2/M transition of mitotic cell cycle; positive regulation of protein phosphorylation; and regulation of gene expression. Located in centrosome; cytosol; and nuclear lumen. Implicated in pre-eclampsia. Biomarker of Alzheimer's disease.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

219736

Swiss Prot

Q6ZVD7

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 520-800 of human STOX1 (NP_689922.3).

Synonyms

C10orf24; STOX1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.