

LYST Rabbit pAb

Catalog No.: A18139

Basic Information

Observed MW

Refer to figures

Calculated MW

429kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse, Rat

Background

This gene encodes a protein that regulates intracellular protein trafficking in endosomes, and may be involved in pigmentation. Mutations in this gene are associated with Chediak-Higashi syndrome, a lysosomal storage disorder. Alternative splicing results in multiple transcript variants, though the full-length nature of some of these variants has not been determined.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

1130

Swiss Prot

Q99698

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human LYST (NP_000072.2).

Synonyms

CHS; CHS1; Mauve; LYST

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.