

MRGPRF Grp168 Rabbit pAb

Catalog No.: A18159

Basic Information

Observed MW

Refer to figures

Calculated MW

38kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Predicted to enable G protein-coupled receptor activity. Predicted to be involved in G protein-coupled receptor signaling pathway. Located in nuclear membrane and plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

116535

Swiss Prot

Q96AM1

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 250-343 of human MRGPRF Grp168 (NP_001091985.1).

Synonyms

RTA; MRGF; GPR140; GPR168; MRGPRF Grp168

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.