

FAM192A Rabbit pAb

Catalog No.: A18172

Basic Information

Observed MW

Refer to figures

Calculated MW

29kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse

Background

Involved in negative regulation of proteasomal protein catabolic process and negative regulation of protein binding activity. Located in nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

80011

Swiss Prot

Q9GZU8

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 120-200 of human FAM192A (NP_079222.1).

Synonyms

CDA10; NIP30; PIP30; CDA018; FAM192A; C16orf94

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.