

TMK1 Rabbit pAb

Catalog No.: A18176

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Arabidopsis thaliana

Background

receptor-like transmembrane kinase I (TMK1)

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

842930

Swiss Prot

Immunogen

Recombinant protein of human TMK1.

Synonyms

F15E12.4; F15E12_4; transmembrane kinase 1; TMK1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.