

SUPT3H Rabbit pAb

Catalog No.: A18179

Basic Information

Observed MW

36kDa

Calculated MW

36kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse

Background

Enables transcription coactivator activity. Involved in histone H3 acetylation and histone deubiquitination. Located in nucleoplasm. Part of SAGA complex and transcription factor TFTC complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

8464

Swiss Prot

O75486

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 50-130 of human SUPT3H (NP_852001.1).

Synonyms

SPT3; SPT3L; SUPT3H

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from A-549 cells, using SUPT3H Rabbit pAb (A18179) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 3min.