

# CFAP69 Rabbit pAb

Catalog No.: A18193

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

106kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse

## Background

Acts upstream of or within sperm axoneme assembly. Located in cytoplasm and sperm midpiece. Implicated in spermatogenic failure 24.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

79846

### Swiss Prot

A5D8W1

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human CFAP69 (NP\_001034795.2).

### Synonyms

FAP69; SPGF24; C7orf63; CFAP69

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.